




We Create Technologies
for Interactive Retail Experiences

Dropper

Dropper is the new way to product sampling with a no-contact, ultra-compact solution that connects directly to a product's original packaging. This innovative device precisely dispenses almost any liquid product—from serums and creams to fragrances and more—ensuring each sample is delivered without waste or mess.

How it works

Dropper interfaces seamlessly with a product's existing container, promoting efficient and precise dispensing. Activation is achieved through a non-contact sensor, making the process hygienic and simple. This method preserves the product's integrity, minimizes waste, and eliminates the mess typically associated with traditional sampling methods, enhancing the customer experience in retail settings.



Why Choose Dropper ?

Scalability

Smoothly manage operations worldwide with our all-in-one maintenance system.

Hygiene

Provides a no-contact, cleaner way to sample cosmetics, and adaptable to various liquid products

Eco-Conscious

Reduces environmental impact by using 5 to 20 times fewer bottles, eliminates theft, and employs a sustainable rental model.

Marketing Data

Collect insights on consumer behavior and engagement to refine marketing strategies and better manage inventory.

SKINCARE FAVORITES

SAMPLE OUR BESTSELLING SKINCARE!



Type of products



Fragrances



Cream



Serum/Oil



Foundation

Installation & maintenance



Worldwide

General information

POWER : 5V & 1,5A

POWER SUPPLY : wire & battery



CB

CERTIFICATION :

- **Standard IEC for CB:** - IEC 60335-1:2010, COR1:2010, COR2:2011, AMD1:2013, COR1:2014, AMD2:2016, COR1:2016
- **Standard for cTUVus:** UL 60335-1 :16 (E

Production capacity : in stocks for small quantities or 3 months of production time.

Type of products adaptable : fragrances, spirits, oils, serums, gels and rich textures (cream, foundation, ...).

Number of tests and device autonomy :

- **Fragrances** : 3000 tests for 100ml (3, 1/3 Oz)
- **Rich moisturizers** : Between 250 and 350 tests for 50ml (1,69 Oz)
- Can also be used with papers or any support

Model types : Horizontal, vertical, multiple.

Option : motorized caps for rich moisturizers.


Managed doses for fragrances and rich textures: several pumps available

- Fragrances : 30 μ L
- Skincare and foundation : 70 μ L, 150 μ L , 250 μ L


Detection height : from 20mm to 180mm

Optional luminous feedback with lumisheet or any led.

Naked Device



Integrated in fixture



Business Model

Machine rental according 3 models:

- Events
- Multiple Events
- Permanent

Optional

Purchase of our software licenses according to 3 models :

- Events
- Multiple Events
- Permanent

**Click here to generate
a cost estimate**



Maintenance

Our maintenance is ensured by GK Cloud, that redefines maintenance with a focus on proactive, reactive, and preventive strategies. Our platform is built on the principle of uninterrupted service, ensuring retail technologies are always operational and effective.

ALL THE FEATURES INCLUDED IN THE MAINTENANCE

- Online Support
- Dropper Technology
- Platform
 - Marketing Data
 - Automated Report
 - Predictive Maintenance Alert
- Mobile app for field
 - Notification center for fields
 - Synchronisation of datas in case there is not Internet access
 - Automated hardware updates on App
- Customer success
 - Proactive customer success management
 - Phone Assistance
 - On site intervention management

Plans

